

Statement

We Getein Biotech, Inc. herein state that the variants Alpha, Beta, Gamma, Delta, and Omicron will not affect diagnosis sensitivity of product **One Step Test for SARS-CoV-2 Antigen (Colloidal Gold)**.

The One Step Test for SARS-CoV-2 Antigen (Colloidal Gold) by our company detects N protein of SARS-CoV-2. While the mutation of the virus found in England (Alpha variant B.1.1.7), South Africa (Beta variant B.1.351), Brazil (Gamma variant P.1), India (Delta variant B.1.617.2), and Multiple countries (Omicron variant B.1.1.529) are not included in the detection sequence of our product and therefore our product is still effective to detect the antigen of novel coronavirus and the mutation will not affect the detection result of our product.

We will keep evaluating the impact of new variants.

